


CUBESCAN™ BIOCON-700 BLADDER SCANNER


de smit medical

Advanced technology



A popular, versatile and fully featured bladder scanner for reliable and accurate results with data transfer to your EPR

Versatility



SIMPLE TO USE

ACCURATE & PROVEN

ROBUST & RELIABLE

EPR CONNECTION

VALUE FOR MONEY


Routinely used



**CHOICE
NOT COMPROMISE**

TECHNOLOGY FIRSTS IN BLADDER SCANNING

Inventors of Pre-Scan - now adopted as industry standard



Advanced software

Bladderpoint™ auto correction to assist all users
Child scan depth marker



Data transfer

Upload patient readings to EPR*
XML/JPG/PDF OR HL7

KEY FEATURES

Simple and universal use

Appropriate for all patient types
Suitable for all levels of user
7" colour touch display

Data security

Patient data encryption

Value for money

Lower cost of purchase and exceptionally low whole life cost
Robust probe design to reduce risk of damage

Reliable

No annual routine calibration is required
Self diagnostic test function

Accurate

Clinically proven accuracy
3D modelling with auto correction and Bladderpoint™ provides accurate readings for all

Extremely robust

View our probe drop test video
Low number of probe returns

TRAINING & TECHNICAL SUPPORT

All our bladder scanners are provided with free, competency based, train-the-trainer user training along with a comprehensive user manual and quick reference guide
A comprehensive virtual training package is available

24 month warranty as standard which includes accidental damage cover on all parts during the first 12 months (T&Cs apply)
Flexible service solutions available along with technical training for in-house servicing and calibration

All technical support is delivered from our UK offices

SPECIFICATION

Accuracy: 0-99ml ± 15 ml, 100-999ml ± 15 %

Range: 0-999ml

Size & weight: H 13cm, W 26cm, D 6.6cm; 1.3kg

Power: lithium-ion battery or mains operated

Battery: to fully charged <6 hours

2400 patient scans from fully charged battery

Accessories: trolley mounted or carry case

Display: 7" colour TFT-LCD

Internal data storage: 40 patient records

External data storage: 1400 patient records

Scan time: 2 seconds with onboard scan guide

Data transfer to EPR: XML/JPG/PDF or HL7

***EPR:** connectivity through CUBEpro PC software via USB

ALSO IN THE FAMILY



de smit medical

Contact Us

0845 345 4226

www.desmitmedical.com

sales@desmitmedical.com

Manufacturer: Mcube Technology